



## HOW RADIOFREQUENCY TREATMENT RELIEVES CHRONIC NECK AND BACK PAIN


1. A medial branch nerve transmits pain signals from the spinal facet joint.



2. A needle is inserted in treatment area. Stimulation is used to locate nerves.



3. Radiofrequency current heats nearby nerve tissue.



4. Pain signals are blocked.



The Cosman Radiofrequency Generators, associated Radiofrequency Lesion Probes and RF Cannula are indicated for use in procedures to create radiofrequency lesions for the treatment of pain or for lesioning nerve tissue for functional neurosurgical procedures. The Cosman RF Injection Electrodes are used for percutaneous nerve blocks with local anesthetic solution or for radiofrequency lesioning. Refer to the Instructions for Use provided with Cosman generators, electrodes and cannulas for potential adverse effects, warnings and precautions prior to using these products. Caution: Federal (U.S.) law restricts this device to sale by or on the order of a physician.

Outside of US Indications for Use: CAUTION: The law restricts these devices to sale by or on the order of a physician. Indications, contraindications, warnings and instructions for use can be found in the product labelling supplied with each device. Information for use only in countries with applicable health authority registrations. Material not intended for use in France.

Products shown for INFORMATION purposes only and may not be approved or for sale in certain countries. Please check availability with your local sales representative or customer service.

Procedures presented are not meant as replacement for the product Directions for Use.

All trademarks are the property of their respective owners.

**Boston Scientific**  
Advancing science for life™

25155 Rye Canyon Loop  
Valencia, CA 91355 USA

Copyright ©2017 by  
Boston Scientific Corporation  
or its affiliates. All rights reserved.

NM-482822-AA\_NOV2017

**Boston Scientific**  
Advancing science for life™

# RADIOFREQUENCY PAIN MANAGEMENT

for Chronic Neck and Back Pain



## PAIN MANAGEMENT WITH RADIOFREQUENCY

You may be one of the many chronic back and neck pain patients who visit pain management specialists every year.

Imagine an outpatient procedure that could treat your neck or back pain, and allow you to return quickly to your normal activities. All this may be possible using the Radiofrequency (RF) technology from Boston Scientific.

## RADIOFREQUENCY TREATMENT

Radiofrequency ablation (RFA) is a safe, proven means of interrupting pain signals, such as those coming from an irritated facet joint in the spine. Radiofrequency current is used to heat up a small volume of nerve tissue, thereby interrupting pain signals from that specific area. Clinical data shows that radiofrequency treatment can effectively provide long-lasting pain relief.<sup>1,2</sup>

Radiofrequency Ablation (RFA) Therapy is a simple and effective treatment for chronic pain originating from the neck, upper back, and lower back. RFA is also known as facet denervation or thermocoagulation. Your doctor can determine if RFA treatment is right for you.

## DURING THE PROCEDURE

A radiofrequency nerve ablation procedure is performed by a doctor in a treatment room setting. Both local anesthesia and a mild sedative may be used to reduce discomfort during the procedure.

During the procedure you will be lying on either your stomach or back, depending on the treatment location. Your doctor will insert a small needle in your neck or back. Using x-ray imaging, the doctor will guide the needle to the treatment location near the spine. An electrode is inserted through the needle to stimulate nearby sensory and motor nerves. During stimulation, your doctor will ask if you are able to feel a tingling sensation. The object of stimulation is to determine the optimal treatment location for pain relief.

Once the needle and electrode placement is verified, a small RF current will travel through the electrode into the surrounding tissue, which is designed to cause the tissue to heat and interrupt pain signals. You should alert your doctor if you experience discomfort at any time during the procedure.

Your doctor may treat more than one location during the same procedure.



This pamphlet is for general education only. All medical questions and concerns should be addressed with your doctor. Your doctor can explain possible outcomes of treatment including risks and side effects.

## AFTER THE PROCEDURE

You may experience some soft tissue discomfort at the needle placement site(s) following the procedure. Like other soft tissue wounds, the discomfort should subside over several days or weeks. If there are no complications, the procedure can usually be done on an outpatient basis.

Your back pain should subside over the next several weeks. Radiofrequency treatment of nerves usually blocks pain signals for a prolonged period of time, up to two years.<sup>2</sup> The procedure can be repeated if the nerves regenerate and your pain returns.

Radiofrequency treatment should not limit your daily activities. Physical restrictions you had prior to the procedure may still remain. As with any medical procedure, however minor, there are certain risks involved. Please ask your doctor for details regarding the potential risks with radiofrequency ablation and what activities are appropriate after radiofrequency treatment.

**Ask your doctor about pain treatment with Boston Scientific radiofrequency equipment.**

1. Gauci CA. Radiofrequency treatment of the lumbar medial branch. Cosman Procedure Technique Series. USA; 2009.

2. Lord SM, et al. Percutaneous radiofrequency for chronic cervical zygapophyseal joint pain The New England Journal of Medicine 1996; 335(23): 1721-1726.